



Recombinant Human Calcipressin-1 (RCAN1)

Product Code	CSB-BP019500HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P53805
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MEEVDLQDLPSATIACHLDPRVFDGLCRAKFESLFRTYDKDITFQYFKSFKRV RINFSPFSAADARLQLHKTEFLGKEMKLYFAQTLHIGSSHLAPPNPKQFLIS PPASPPVGWKQVEDATPVINYDLLYAISKLGPEKEYELHAATDTTPSVVVHVC ESDQEKEEEEEEMERMRRPKPKIIQTRRPEYTPIHLS
Source	Baculovirus
Target Names	RCAN1
Protein Names	Recommended name: Calcipressin-1 Alternative name(s): Adapt78 Down syndrome critical region protein 1 Myocyte-enriched calcineurin-interacting protein 1 Short name= MCIP1 Regulator of calcineurin 1
Expression Region	1-197
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length of isoform 2
Target Details	This protein interacts with calcineurin A and inhibits calcineurin-dependent signaling pathways, possibly affecting central nervous system development. This gene is located in the minimal candidate region for the Down syndrome phenotype, and is overexpressed in the brain of Down syndrome fetuses. Chronic overexpression of this gene may lead to neurofibrillary tangles such as those associated with Alzheimer disease. Three transcript variants encoding three different isoforms have been found for this gene.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.