



LDLRAP1 Antibody

Product Code	CSB-PA012850GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5SW96
Immunogen	Human LDLRAP1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	low density lipoprotein receptor adaptor protein 1;LDLRAP1;ARH;ARH1;ARH2;DKFZp586D0624;FHCB1;FHCB2;MGC34705 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	LDLRAP1
Target Details	This protein is a cytosolic protein which contains a phosphotyrosine binding (PTD) domain. The PTD domain has been found to interact with the cytoplasmic tail of the LDL receptor. Mutations in this gene lead to LDL receptor malfunction and cause the disorder autosomal recessive hypercholesterolaemia.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.