



LILRA2 Antibody

Product Code	CSB-PA012934GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N149
Immunogen	Human LILRA2
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	leukocyte immunoglobulin-like receptor, subfamily A (with TM domain), member 2; LILRA2; CD85H; ILT1; LIR-7; LIR7 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	LILRA2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.