



LMBR1L Antibody

Product Code	CSB-PA012996GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6UX01
Immunogen	Human LMBR1L
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	489aa,55kd;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	LMBR1L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.