



LRG1 Antibody

Product Code	CSB-PA013082GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02750
Immunogen	Human LRG1
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	leucine-rich alpha-2-glycoprotein 1;LRG1;HMFT1766;LRG ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	LRG1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.