



LRRC8C Antibody

Product Code	CSB-PA013170GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TDW0
Immunogen	Human LRRC8C
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	leucine rich repeat containing 8 family, member C;LRRC8C;AD158;DKFZp586J1119;FAD158;MGC138551 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	LRRC8C
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.