



LSM2 Antibody

Product Code	CSB-PA013205GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y333
Immunogen	Human LSM2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	LSM2 homolog, U6 small nuclear RNA associated (S. cerevisiae);LSM2;C6orf28;G7b;YBL026W;snRNP ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	LSM2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.