



MAGEH1 Antibody

Product Code	CSB-PA013358GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H213
Immunogen	Human MAGEH1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	melanoma antigen family H, 1;MAGEH1;APR-1;APR1;MAGEH ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MAGEH1
Target Details	This gene is thought to be involved in apoptosis. Multiple polyadenylation sites have been found for this gene.