



MAP1S Antibody

Product Code	CSB-PA013407GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q66K74
Immunogen	Human MAP1S
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB,IHC,IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	microtubule-associated protein 1S;MAP1S;BPY2IP1;C19orf5;FLJ10669;MAP8;MGC133087;VCY2IP-1;VCY2IP 1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MAP1S
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.