



# MAPKAP1 Antibody

<b>Product Code</b>	CSB-PA013472GA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9BPZ7
<b>Immunogen</b>	Human MAPKAP1
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	ELISA,WB,IHC,IF
<b>Storage Buffer</b>	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
<b>Purification Method</b>	Antigen Affinity purified
<b>Isotype</b>	IgG
<b>Alias</b>	mitogen-activated protein kinase associated protein 1;MAPKAP1;JC310;MGC2745;MIP1;SIN1;SIN1b;SIN1g ;
<b>Product Type</b>	Purified Rabbit Anti human PolyClonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	MAPKAP1
<b>Target Details</b>	This gene encodes a protein that is highly similar to the yeast SIN1 protein, a stress-activated protein kinase. Alternatively spliced transcript variants encoding distinct isoforms have been described. Alternate polyadenylation sites as well as alternate 3 UTRs have been identified for transcripts of this gene.
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.