



# MARCKSL1 Antibody

<b>Product Code</b>	CSB-PA013494GA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P49006
<b>Immunogen</b>	Human Marcksl1
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse
<b>Tested Applications</b>	ELISA,WB
<b>Storage Buffer</b>	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Alias</b>	MARCKS-like 1;Marcksl1;AL022768;AW215397;AW536807;D4Bc1;F52;MacMARCKS;Macs-2 ;Macs-3;Macs2;Macs3;Mlp;Mrp ;
<b>Product Type</b>	Purified Rabbit Anti mouse PolyClonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	MARCKSL1
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.