



MARS Antibody

Product Code	CSB-PA013500GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P56192
Immunogen	Human MARS
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	methionyl-tRNA synthetase;MARS;FLJ35667;METRS;MTRNS ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MARS
Target Details	Aminoacyl-tRNA synthetases are a class of enzymes that charge tRNAs with their cognate amino acids. This protein belongs to the class I family of tRNA synthetases.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.