



MAVS Antibody

Product Code	CSB-PA013526GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z434
Immunogen	Human MAVS
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB,IF
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	mitochondrial antiviral signaling protein;MAVS;CARDIF;DKFZp547C224;DKFZp666M015;FLJ27482;FLJ35386;FLJ38051;FLJ41962;IPS-1;IPS1;KIAA1271;MGC3260;VISA ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MAVS
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.