



# Recombinant Clostridium kluveri Elongation factor Ts (tsf)

<b>Product Code</b>	CSB-YP025124CWX
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	A5N829
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Clostridium kluveri (strain ATCC 8527 / DSM 555 / NCIMB 10680)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MITAQMVKEL RERTGAGMMD CKKALNEAGG DSEKAIEILR EKGLAAAANK SGRIASEGLV KTYISEDGKV ASIVEVNCET DFVAVNADFV NFVDNLAKQI SLSESTTVEE LSEEKYISDD SKTVSETLVN LISKLGENMA IRRFERLAVS KGLIESYIHG GGRIGVLVKL ECEKESEILK EVAKDVAMQV AATNPLFLSK DTVDSATLDK EKEIFKVQAL NEGKPEKIAE KIVIGRVQKY YKENCLIEQL WVKDSLTDID KYLKSKSKEV GAPIKISNFI RFEKGEGIEK KEEDFAEEVR KQIEGK
<b>Source</b>	Yeast
<b>Target Names</b>	tsf
<b>Protein Names</b>	Recommended name: Elongation factor Ts Short name= EF-Ts
<b>Expression Region</b>	1-306
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.