



MED31 Antibody

Product Code	CSB-PA013664GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y3C7
Immunogen	Human MED31
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	mediator complex subunit 31;MED31;3110004H13Rik;CGI-125;FLJ27436;FLJ36714;Soh1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MED31
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.