



Recombinant Legionella pneumophila Succinyl-CoA ligase [ADP-forming] subunit beta (sucC)

Product Code	CSB-EP022917LLJ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A5IH21
Product Type	Recombinant Protein
Immunogen Species	Legionella pneumophila (strain Corby)
Purity	≥85% (SDS-PAGE)
Sequence	MNLHEYQAKQ LFASYGLPVP RGEVAYNVED ALLVASQLST SRWVVKAQVH AGGRGKAGGV KLVSSKDELA AVAKSMLGTR LVTYQTDARG QPVNAILVEE TCEIDKELYL GAVVDRSTRR VVIMASTEGG VEIEKVAHET PEKIFKVVVD PLVGVMPFQC RETAFKLGLK DDQIKQFTHL MMGLGKMFVD CDLSLLEINP LVITKSGQLI CLDGKINIDG NALFRQPKLK NMRDVSQEDD RENRASDWEL NYIPLDGTIG CMVNGAGLAM ATMDVIKLHG GEPANFLDVG GGATKERVSE ALKIIVSDEK VKGILVNIFG GIVRCDLIAD GILAAVKEVD VKIPVVVRLE GNNAQLGAEI LNKSNNLVIA ATSLTDAACK IVAAVSE
Source	E.coli
Target Names	sucC
Protein Names	Recommended name: Succinyl-CoA ligase [ADP-forming] subunit beta EC= 6.2.1.5 Alternative name(s): Succinyl-CoA synthetase subunit beta Short name= SCS-beta
Expression Region	1-387
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.