



MGEA5 Antibody

Product Code	CSB-PA013786GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60502
Immunogen	Human MGEA5
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC,IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	meningioma expressed antigen 5 (hyaluronidase);MGEA5;FLJ11229;FLJ23355;KIAA0679;MEA5;NCOAT;OGA ;O-GlcNAcase
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MGEA5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.