



MICAL2 Antibody

Product Code	CSB-PA013808GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O94851
Immunogen	Human MICAL2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,IHC,IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	microtubule associated monooxygenase, calponin and LIM domain containing 2;MICAL2;DKFZp686H03148;DKFZp686H2469;KIAA0750;MICAL2PV1;MICAL2PV2 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MICAL2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.