



MIER1 Antibody

Product Code	CSB-PA013823GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N108
Immunogen	Human MIER1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	mesoderm induction early response 1 homolog (Xenopus laevis);MIER1;DKFZp781G0451;ER1;KIAA1610;MGC131940;MGC150640;MGC150641;MI-ER1;RP5-944N15.1;hMI-ER1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MIER1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.