



MLN Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA014638GA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P12872 |
| Immunogen | Human MLN |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA, WB |
| Storage Buffer | PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles. |
| Purification Method | Antigen Affinity purified |
| Isotype | IgG |
| Alias | motilin;MLN;MGC138519 ; |
| Product Type | Purified Rabbit Anti Human PolyClonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | MLN |
| Target Details | This gene encodes a small peptide hormone that is secreted by cells of the small intestine to regulate gastrointestinal contractions and motility. Proteolytic processing of the secreted protein produces the mature peptide and a byproduct referred to as motilin-associated peptide (MAP). Two transcript variants encoding different preproprotein isoforms but the same mature peptide have been found for this gene. |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |