



MMS19 Antibody

Product Code	CSB-PA014682GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96T76
Immunogen	Human MMS19
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	MMS19 nucleotide excision repair homolog (S. cerevisiae);MMS19;FLJ34167;FLJ95146;MET18;MGC99604;MMS19L;hMMS19 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MMS19
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.