



MND1 Antibody

Product Code	CSB-PA014687GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BWT6
Immunogen	Human MND1
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purifie
Isotype	IgG
Alias	meiotic nuclear divisions 1 homolog (S. cerevisiae);MND1;GAJ ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MND1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.