



MOB1B Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA014695GA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q7L9L4 |
| Immunogen | Human MOBKL1A |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse,Rat |
| Tested Applications | ELISA,WB,IHC |
| Storage Buffer | PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles. |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Alias | MOB1, Mps One Binder kinase activator-like 1A (yeast);MOBKL1A;MATS2;MGC33910;MOB4A;Mob1B ; |
| Product Type | Purified Rabbit Anti human PolyClonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | MOB1B |
| Target Details | This protein is similar to the yeast Mob1 protein. Yeast Mob1 binds Mps1p, a protein kinase essential for spindle pole body duplication and mitotic checkpoint regulation. |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |