



MRI1 Antibody

Product Code	CSB-PA014799GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BV20
Immunogen	Human MRI1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	methylthioribose-1-phosphate isomerase homolog (S. cerevisiae);MRI1;MGC3207;MTNA;Ypr118w ;M1Pi; MRDI
Product Type	Purified Rabbit Anti Human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MRI1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.