



MRPL36 Antibody

Product Code	CSB-PA014843GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P0J6
Immunogen	Human MRPL36
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Form	Liquid
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	Antigen Affinity purified
Isotype	IgG
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MRPL36
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.