



MS4A6A Antibody

Product Code	CSB-PA015018GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H2W1
Immunogen	Human MS4A6A
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	membrane-spanning 4-domains, subfamily A, member 6A;MS4A6A;4SPAN3;4SPAN3.2;CD20L3;CDA01;MGC131944;MGC22650;MS4A6;MST090;MSTP090 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MS4A6A
Target Details	This gene encodes a member of the membrane-spanning 4A gene family. Members of this nascent protein family are characterized by common structural features and similar intron/exon splice boundaries and display unique expression patterns among hematopoietic cells and nonlymphoid tissues. The gene encoding this protein is localized to 11q12.1, among a cluster of family members. Alternative splicing of this gene results in several transcript variants.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.