



MSRB2 Antibody

Product Code	CSB-PA015052GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y3D2
Immunogen	Human MSRB2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	methionine sulfoxide reductase B2;MSRB2;CBS-1;CBS1;CGI-131;MGC26104;MSRB;PILB ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MSRB2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.