



MTIF3 Antibody

Product Code	CSB-PA015166GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H2K0
Immunogen	Human MTIF3
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	mitochondrial translational initiation factor 3; MTIF3; IF-3mt; IF3(mt) ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MTIF3
Target Details	This gene encodes a translation initiation factor that is involved in mitochondrial protein synthesis. Polymorphism in this gene is associated with the onset of Parkinson s disease. Alternate splicing results in multiple transcript variants. A pseudogene of this gene is found on chromosome 5.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.