



# FAN1 Antibody

<b>Product Code</b>	CSB-PA015179GA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9Y2M0
<b>Immunogen</b>	Human MTMR15
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, WB
<b>Storage Buffer</b>	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
<b>Purification Method</b>	Antigen Affinity purified
<b>Isotype</b>	IgG
<b>Alias</b>	myotubularin related protein 15;MTMR15;DKFZp451H236;DKFZp686K16147;KIAA1018 ;
<b>Product Type</b>	Purified Rabbit Anti human PolyClonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	FAN1
<b>Target Details</b>	This gene encodes a member of the myotubularin-related class 1 cysteine-based protein tyrosine phosphatases. The encoded protein may be catalytically inactive. Alternatively spliced transcript variants have been described.
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.