



# MUM1 Antibody

<b>Product Code</b>	CSB-PA015236GA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q2TAK8
<b>Immunogen</b>	Human MUM1
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, WB, IHC
<b>Storage Buffer</b>	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Alias</b>	melanoma associated antigen (mutated) 1;MUM1;FLJ14868;FLJ22283;HSPC211;MGC131891;MGC163315;MUM-1 ;
<b>Product Type</b>	Purified Rabbit Anti human PolyClonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	MUM1
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.