



MYBBP1A Antibody

Product Code	CSB-PA015262GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BQG0
Immunogen	Human MYBBP1A
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC, IF
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	MYB binding protein (P160) 1a; MYBBP1A; FLJ37886; P160; PAP2 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MYBBP1A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.