



Recombinant Bovine Growth/differentiation factor 6 (GDF6)

Product Code	CSB-MP009350BO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P55106
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	≥85% (SDS-PAGE)
Sequence	TAFA SRHGKRHGKK SRLRCSKKPL HVNFKELGWD DWIIPLEYE AYHCEGVCDF PLRSHLEPTN HAIQTLMNS MDPGSTPPSC CVPTKLTPIS ILYIDAGNNV VYNEYEEMVV ESCGCR
Source	Mammalian cell
Target Names	GDF6
Protein Names	Recommended name: Growth/differentiation factor 6 Short name= GDF-6 Alternative name(s): Cartilage-derived morphogenetic protein 2 Short name= CDMP-2
Expression Region	317-436
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene encodes a member of the bone morphogenetic protein (BMP) family and the TGF-beta superfamily of secreted signaling molecules. It is required for normal formation of some bones and joints in the limbs, skull, and axial skeleton. Mutations in this gene result in colobomata, which are congenital abnormalities in ocular development, and in Klippel-Feil syndrome (KFS), which is a congenital disorder of spinal segmentation.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.