



MYOCD Antibody

Product Code	CSB-PA015357GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IZQ8
Immunogen	Human Myocd
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	myocardin;Myocd;BSAC2A;Srfcp ;
Product Type	Purified Rabbit Anti Mus PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MYOCD