



N6AMT2 Antibody

Product Code	CSB-PA015395GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WVE0
Immunogen	Human N6AMT2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	N-6 adenine-specific DNA methyltransferase 2 (putative);N6AMT2;ESP13 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	EEF1AKMT1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.