



Recombinant *Flavobacterium johnsoniae* Lipoyl synthase (lipA)

Product Code	CSB-EP012927FDT
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A5FNF9
Product Type	Recombinant Protein
Immunogen Species	<i>Flavobacterium johnsoniae</i> (strain ATCC 17061 / DSM 2064 / UW101) (<i>Cytophaga johnsonae</i>)
Purity	>85% (SDS-PAGE)
Sequence	METVIENIPP AKPKWLKVKL PIGQKYTELRLGLVDKYSLNT ICTSGSCPNNM GECWGEGTAT FMILGNVCTR SCGFCGVKTG RPETVDWDEP EKVARSIKIM NIKHAVITSV DRDDLKDGGS IIWMETVRAI RRMNPNTTLE TLIPDFQGIE RNIDRIVEAN PEVVSHNMET VRLRTREVRI QAKYDRSLEV LRYLKEKGIN RTKSGIMLGL GETEEEVYQT MRDLRDANVD VVTIGQYLQP TTKHLPVKEF ITPELFAKYE KYGIELGFRH VESGPLVRSS YKAQKHIL
Source	<i>E. coli</i>
Target Names	lipA
Protein Names	Recommended name: Lipoyl synthase EC= 2.8.1.8 Alternative name(s): Lip-syn Short name= LS Lipoate synthase Lipoic acid synthase Sulfur insertion protein LipA
Expression Region	1-288
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.