



NAT10 Antibody

Product Code	CSB-PA015465GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H0A0
Immunogen	Human NAT10
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IF
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	N-acetyltransferase 10 (GCN5-related);NAT10;ALP;DKFZp434C116;FLJ10774;FLJ12179;FLJ23850;KIAA1709 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NAT10
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.