



NCKAP1 Antibody

Product Code	CSB-PA015532GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2A7
Immunogen	Human NCKAP1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	NCK-associated protein 1;NCKAP1;FLJ11291;HEM2;KIAA0587;MGC8981;NAP1;NAP125 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NCKAP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.