



NDUFAF1 Antibody

Product Code	CSB-PA015637GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y375
Immunogen	Human NDUFAF1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	NADH dehydrogenase (ubiquinone) 1 alpha subcomplex, assembly factor 1;NDUFAF1;CGI-65;CGI65;CIA30 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NDUFAF1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.