



NDUFB11 Antibody

Product Code	CSB-PA015643GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NX14
Immunogen	Human NDUFB11
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	NADH dehydrogenase (ubiquinone) 1 beta subcomplex, 11, 17.3kDa; NDUFB11; ESSS; FLJ20494; MGC111182; NP17.3; Np15; P17.3 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NDUFB11
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.