



NDUFB9 Antibody

Product Code	CSB-PA015657GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y6M9
Immunogen	Human NDUFB9
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	NADH dehydrogenase (ubiquinone) 1 beta subcomplex, 9, 22kDa;NDUFB9;B22;DKFZp566O173;FLJ22885;LYRM3;UQOR22 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NDUFB9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.