



NEDD4L Antibody

Product Code	CSB-PA015682GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96PU5
Immunogen	Human NEDD4L
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	neural precursor cell expressed, developmentally down-regulated 4-like;NEDD4L;FLJ33870;KIAA0439;NEDD4-2;RSP5;hNedd4-2 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NEDD4L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.