



NEK9 Antibody

Product Code	CSB-PA015707GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TD19
Immunogen	Human NEK9
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC,IF
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	NIMA (never in mitosis gene a)- related kinase 9;NEK9;DKFZp434D0935;MGC138306;MGC16714;NERCC;NERCC1;Nek8 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NEK9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.