



NFIA Antibody

Product Code	CSB-PA015754GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12857
Immunogen	Human NFIA
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA, WB, IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	nuclear factor I/A; NFIA; DKFZp434L0422; DKFZp686J23256; FLJ39164; KIAA1439; NFI-L ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NFIA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.