



NFKB2 Antibody

Product Code	CSB-PA015760GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q00653
Immunogen	Human NFKB2
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	nuclear factor of kappa light polypeptide gene enhancer in B-cells 2 (p49/p100); NFKB2; LYT-10; LYT10 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NFKB2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.