



NIF3L1 Antibody

Product Code	CSB-PA015806GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9GZT8
Immunogen	Human NIF3L1
Raised In	Rabbit
Species Reactivity	Human,Rat
Tested Applications	ELISA,WB
Form	Liquid
Storage Buffer	PBS with 0.02% sodium azide and 50% glycerol, pH 7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NIF3L1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.