



NIPSNAP3A Antibody

Product Code	CSB-PA015819GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UFN0
Immunogen	Human NIPSNAP3A
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	nipsnap homolog 3A (C. elegans);NIPSNAP3A;DKFZp564D177;FLJ13953;HSPC299;MGC14553;NIPSNAP4;TASSC ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NIPSNAP3A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.