



NLRP2 Antibody

Product Code	CSB-PA015870GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NX02
Immunogen	Human NLRP2
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	NLR family, pyrin domain containing 2;NLRP2;CLR19.9;FLJ20510;NALP2;NBS1;PAN1;PYPAF2 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NLRP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.