



NME5 Antibody

Product Code	CSB-PA015890GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P56597
Immunogen	Human NME5
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	non-metastatic cells 5, protein expressed in (nucleoside-diphosphate kinase);NME5;NM23-H5;NM23H5 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NME5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.