



NMUR1 Antibody

Product Code	CSB-PA015903GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HB89
Immunogen	Human NMUR1
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	neuromedin U receptor 1; NMUR1; (FM-3); FM-3; FM3; GPC-R; GPR66; NMU1R ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NMUR1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.