



NPL Antibody

Product Code	CSB-PA015994GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BXD5
Immunogen	Human NPL
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	N-acetylneuraminase pyruvate lyase (dihydrodipicolinate synthase);NPL;C1orf13;MGC149582;MGC61869;NPL1;c112 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NPL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.